


Simply faster. Simply more reliable.




T-REF Chamber


Reference Detector
for Small Fields

T-REF Chamber

The New Reference Detector for Small Field Relative Dosimetry



Easy to mount on edge of water tank



Large diameter to avoid positioning problems

When measurement time and accuracy matter, T-REF is the answer.

With the new T-REF chamber, small field PDD and profile measurements can be performed much faster and more reliably.

As a thin, plane-parallel transmission chamber that shows a high nominal response and causes no measurable perturbation of the beam, the new T-REF chamber has been designed to avoid positioning problems encountered when using a reference detector for small field dosimetry.


Mounted on the edge of the water tank rather than on a tray connected to the linac head, it ensures precise alignment of the reference detector at any fixed position above 20 cm detector surface distance without perturbing the radiation field or shadowing the field detector.

*Further Information: D. Ceska:
Reference detector for small fields – the T-REF
chamber, Medical Physics International Journal,
published on www.mpjournal.org April 2016*

Highlights

- ▶ Thin-window transmission chamber with large collecting volume (10.5 cm³)
- ▶ Very low total area density of 206 mg/cm²
- ▶ No measurable perturbation of the beam
- ▶ Plug and measure – no repositioning required
- ▶ High and very stable signal
- ▶ No influence by vibrations
- ▶ Excellent signal-to-noise ratio (SNR)
- ▶ No contact to linac head – no accessory tray needed, no problems with temperature drifts
- ▶ Fast and easy to mount
- ▶ Covers entire range of use for small field sizes

Provides reference signal while the beam transmits through the chamber.



WATER SURFACE

SMALL FIELD

FIELD DETECTOR

Very stable signal




Fig. 1: Reference signal of a 4 cm x 4 cm field, 6 MV Varian Truebeam, measured with a Semiflex 0.125 cm³ (T31010) and the T-REF chamber (T34091).

There are different techniques to measure the reference signal. One is to place a Semiflex 0.125 cm³ in the corner of the field. When comparing this technique with the transmission measurement of the reference signal by the

T-REF chamber, the differences can be clearly seen: The signal of the T-REF chamber is more stable (see Fig. 1). Furthermore vibrations would not influence the reference signal.

Perturbation free




Fig. 2: First few millimeters of percentage depth dose curves of a 2 cm x 2 cm field, 6 MV Elekta Synergy, field detector Diode E (T60017), reference detector T-REF chamber (T34091), different distances from 22 cm to 35 cm.

Above 20 cm detector surface distance (DSD) no influence like partial build up effect, scattered radiation or beam hardening is observable in the data (see Fig. 2).

As a result, the T-REF chamber can be placed at any position above 20 cm DSD.

Specifications

T-REF Chamber Type 34091

| | |
|--------------------------|---|
| Type of product | vented plane-parallel ionization chamber |
| Application | relative dosimetry in high-energy photon beams |
| Nominal sensitive volume | 10.5 cm ³ |
| Design | waterproof, vented, guarded, perturbation-free |
| Reference point | inside of entrance window, center |
| Direction of incidence | perpendicular to the entrance window, see label "Focus" |
| Nominal response | 325 nC/Gy (at ⁶⁰ Co free in air) |
| Chamber voltage | 400 V nominal ± 500 V maximal |
| Polarity effect | ≤ ± 1 % |
| Leakage current | ≤ ± 100 fA |
| Cable leakage | ≤ 1 pC/(Gy·cm) |

Materials and measures

| | |
|-----------------------------------|----------------------------|
| Total area density | 206 mg/cm ² |
| Water-equivalent window thickness | 2.06 mm for photons |
| Transmission | > 99 % for energies ≥ 6 MV |
| Dimension of sensitive volume | radius 40.8 mm, depth 2 mm |

Ion collection efficiency at nominal voltage

| | |
|-------------------------|---------|
| Ion collection time | 67 μs |
| Max. dose rate for | |
| ≥ 99.5 % saturation | 21 Gy/s |
| ≥ 99.0 % saturation | 42 Gy/s |
| Max. dose per pulse for | |
| ≥ 99.5 % saturation | 0.9 mGy |
| ≥ 99.0 % saturation | 1.8 mGy |

Useful ranges

| | |
|--|--|
| Chamber voltage | ± (300 ... 500) V |
| Radiation quality | ⁶⁰ Co ... 25 MV photons |
| Max. field size in 20 cm distance to water surface | (5 x 5) cm ² |
| Temperature | (10 ... 40) °C (50 ... 104) °F |
| Humidity | (10 ... 80) %, max 20 g/m ³ |
| Air pressure | (700 ... 1060) hPa |

Ordering Information

| | |
|---------|--|
| TN34091 | T-REF chamber, connecting system BNT, including holder |
| TW34091 | T-REF chamber, connecting system TNC, including holder |
| TM34091 | T-REF chamber, connecting system M, including holder |

Dosimetry Pioneers since 1922.

It all started with a brilliant invention – the revolutionary Hammer dosimeter in 1922. Ingenuity coupled with German engineering know-how shaped the company's history and led to innovative dosimetry solutions that later became an industry standard. Over the years, PTW has maintained its pioneering spirit, growing into a global market leader of dosimetry applications well known for its outstanding quality and precision. Today, PTW dosimetry is one of the first choices for healthcare professionals in radiation therapy, diagnostic radiology, nuclear medicine and health physics.

For more information on PTW dosimetry products, visit www.ptw.de or contact your local PTW representative:

Headquarters

PTW-Freiburg
Physikalisch-Technische
Werkstätten
Dr. Pychlau GmbH
Tel. +49 761 49055-0
info@ptw.de
www.ptw.de

PTW Dosimetría Iberia S. L.
Tel. +34 96 346 2854
info@ptwdi.es
www.ptwdi.es

PTW-Asia Pacific Ltd.
Tel. +852 2369 9234
info@ptw-asiapacific.com
www.ptw-asiapacific.com

PTW-UK Ltd.
Tel. +44 1476 577503
sales@ptw-uk.com
www.ptw-uk.com

PTW-New York Corporation
Tel. (1-516) 827 3181
ptw@ptwny.com
www.ptwny.com

PTW-Beijing Ltd
Tel. +86 10 6443 0746
info@ptw-beijing.com
www.ptw-beijing.com

PTW-France SARL
Tel. +33 1 64 49 98 58
info@ptw-france.com
www.ptw-france.com

PTW-Latin America
Tel. +55 21 2178 2188
info@ptw.com.br
www.ptw.com.br

PTW Dosimetry India Pvt. Ltd.
Tel. +91 44 42079999
info@ptw-india.in
www.ptw-india.in