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## RQT – New CT Radiation Quality

The international standard IEC 61267 “Medical diagnostic X-ray equipment - Radiation conditions for use in the determination of characteristics“ defines the requirements for different diagnostic radiation qualities. The radiation qualities RQT 8, RQT 9 and RQT 10 were defined for CT applications [1]:

Standard RADIATION QUALITY Characterization	X-RAY TUBE VOLTAGE kV	Nominal ADDED FILTER for RQT thickness in copper mm	First HALF-VALUE LAYER in thickness of aluminium mm
RQT 8	100	0,2	6,9
RQT 9	120	0,25	8,4
RQT 10	150	0,3	10,1

In the past PTW calibrated CT ionization chambers at the diagnostic radiation qualities RQR and RQA according to the requirements of the product standard DIN EN 61674 “Medical electrical equipment - Dosimeters with ionization chambers and/or semi-conductor detectors as used in X-ray diagnostic imaging” [2]. Now, the currently available draft of the new product standard IEC 61674 [3] refers to the IEC 61267 and requires RQT radiation qualities for CT applications. Moreover, the revised draft does no longer provide two different calibration factors for attenuated and unattenuated measurements for CT applications. The German national metrology institute *Physikalisch Technische Bundesanstalt* (PTB) in Braunschweig already provides RQT radiation qualities according to the requirements of the revised edition.

To ensure traceability to the IEC and PTB, PTW-Freiburg changed the radiation quality for CT calibration to RQT (reference condition RQT 9) accordingly. Please note the following remarks:

- ▶ The CT ion chamber can be used as usual for RQR 6 to RQR 10 and RQA 3 to RQA 10.
- ▶ The CT ion chamber can be used for RQT 8 to RQT 10 measurements, additionally.
- ▶ The calibration factor for RQT (reference point RQT 9) may differ for approx.  $\pm 4\%$  from the previous RQR calibration factor and by approx.  $\pm 3\%$  from the RQA calibration factor. Therefore the measuring results might deviate from previous measurement results accordingly.
- ▶ The calibration change does not affect the required measuring accuracy acc. to IEC 61674.
- ▶ DIADOS E contains the same calibration factor for attenuated and unattenuated measurements.

Thus pressing the ALT+APPL. keys is no longer required.

- [1] IEC 61267 (2005): Medical diagnostic X-ray equipment - Radiation conditions for use in the determination of characteristics  
 [2] DIN EN 61674 (1997): Medical electrical equipment - Dosimeters with ionization chambers and/or semi-conductor detectors as used in X-ray diagnostic imaging  
 [3] Entwurf/Draft IEC 61674 (2009): Medical electrical equipment - Dosimeters with ionization chambers and/or semi-conductor detectors as used in X-ray diagnostic imaging