


| Treatment Delivery System | Recommended PTW Detectors | MP3 | MP3-M | MP3-T | MP3-XS SRS | MP3-XS IORT | MP3-P | MP2 | MP1 |
|---|--|--|---|---|---|---|---|---|--------|
| | | ¹ L981160 BNT L981161 TNC L981162 M | L981166 BNT L981167 TNC L981168 M | L981199 BNT L981183 TNC L981182 M | L981172 BNT L981173 TNC L981174 M | L981383 BNT L981384 TNC L981385 M | L981185 BNT L981186 TNC L981233 M | L981175 BNT L981176 TNC L981177 M | T41019 |
| Standard LINACS (up to 40 x 40 cm ²) | Semiflex 0.125 cm ³ ; Adv. Markus 0.02 cm ³ ; Roos 0.35 cm ³ ; Diodes P, E, SRS ² | ■ | ■ | - | - | - | - | ■ | ■ |
| TomoTherapy[®] Hi-Art[®] | SFD Mammo 6 cm ³ ; PinPoint 0.03 cm ³ | - | - | ■ | - | - | - | - | - |
| CyberKnife[®] | Semiflex 0.125 cm ³ ; Farmer 0.6 cm ³ ; Diodes E, SRS PinPoint 3D 0.016 cm ³ | - | ■ | - | ■ | - | - | - | - |
| NOVALIS[®] | Semiflex 0.125 cm ³ ; Diodes E, SRS ² | - | ■ | - | - | - | - | - | - |
| BEAM MODULATOR[™] | Semiflex 0.125 cm ³ ; Diodes E, SRS ² | - | ■ | - | - | - | - | - | - |
| VERO SBRT | Semiflex 0.125 cm ³ ; Diodes E, SRS ² | - | ■ | - | - | - | - | - | - |
| IORT Systems | Semiflex 0.125 cm ³ ; Roos 0.35 cm ³ | - | - | - | - | ■ | - | - | - |
| Particle Therapy Systems | Thin-window Bragg Peak, Bragg Peak, 10.5 cm ³ | - | - | - | - | - | ■ | - | - |
| Cobalt Units | Semiflex 0.125 cm ³ | - | ■ | - | - | - | - | ■ | - |

■ recommended or suitable - not available and/or not recommended


¹Part Nos for 230 V, Part Nos for 115 V: L981163, 64, 65 (BNT, TNC, M); ²small photon fields with max. 6 MV



| Key Components | MP3 | MP3-M | MP3-T | MP3-XS SRS | MP3-XS IORT | MP3-P | MP2 | MP1 |
|--|--|---|---|---|---|---|---|--------------|
| | ¹ L981160 BNT L981161 TNC L981162 M | L981166 BNT L981167 TNC L981168 M | L981199 BNT L981183 TNC L981182 M | L981172 BNT L981173 TNC L981174 M | L981383 BNT L981384 TNC L981385 M | L981185 BNT L981186 TNC L981233 M | L981175 BNT L981176 TNC L981177 M | T41019 |
| Scanning Range (cm ³) | 60 x 50 x 40.8 | 50 x 50 x 40.8 | 52 x 15 x 22 | 20 x 20 x 30 | 20 x 20 x 30 | 35 x 38 x 25 | 50 x 40 x 40.7 | 32 x 32 x 37 |
| Motorized Water Tank | 3D | 3D | 3D | 3D | 3D | 3D | 2D | 1D |
| Positioning & Levelling Device | ■ | ■ | - | - | - | - | ■ | - |
| Leg Levellers | - | - | ■ | ■ | ■ | ■ | - | ■ |
| TBA Control Unit & Pendant | ■ | ■ | ■ | ■ | ■ | ■ | ■ | □ |
| TRUFIX [®] Detector Positioning | ■ | ■ | - | □ | □ | - | □ | - |
| TANDEM Electrometer | ■ | ■ | ■ | ■ | - | - | ■ | □ |
| TANDEM ^{xdr} Electrometer | - | - | - | - | ■ | ■ | - | - |
| Field Detector | ■ | ■ | ■ | ■ | ■ | □ | ■ | □ |
| Reference Detector | ■ | ■ | ■ | ■ | ■ | □ | ■ | □ |
| Lifting Carriage | ■ | ■ | - | - | - | - | □ | - |
| Built-in Water Reservoir | - | ■ | - | - | - | - | □ | - |
| Separate Water Reservoir | ■ | - | - | - | - | ■ | - | - |
| MEPHYSTO Software | ■ | ■ | ■ | ■ | ■ | ■ | ■ | □ |

■ included □ optional - not available and/or not recommended


¹Part Nos for 230 V, Part Nos for 115 V: L981163, 64, 65 (BNT, TNC, M)



| Options & Accessories | MP3 | MP3-M | MP3-T | MP3-XS SRS | MP3-XS IORT | MP3-P | MP2 | MP1 |
|--------------------------------|--|---|---|---|---|---|---|--------|
| | ¹ L981160 BNT L981161 TNC L981162 M | L981166 BNT L981167 TNC L981168 M | L981199 BNT L981183 TNC L981182 M | L981172 BNT L981173 TNC L981174 M | L981383 BNT L981384 TNC L981385 M | L981185 BNT L981186 TNC L981233 M | L981175 BNT L981176 TNC L981177 M | T41019 |
| Options | | | | | | | | |
| TRUFIX® Detector Positioning | ■ | ■ | - | □ | □ | - | □ | - |
| LA48 Linear Ion Chamber Array | □ | □ | - | - | - | - | □ | - |
| MLC QA Software for LA48 | □ | □ | - | - | - | - | □ | - |
| PLEXITOM® Verification Phantom | □ | □ | - | □ | - | - | □ | - |
| TPR Measurement | □ | □ | - | - | - | - | □ | - |
| Absolute Dosimetry | ■ | ■ | □ | ■ | ■ | □ | ■ | □ |
| Film Dosimetry | □ | □ | □ | □ | □ | □ | □ | □ |
| Data Transfer to TPS | ■ | ■ | □ | □ | □ | □ | ■ | - |
| MP1 Water Tank (TG-51) | □ | □ | - | - | - | - | - | ■ |
| Thin Entrance Window | □ | □ | - | - | - | - | - | - |
| 3D Detector Block | - | - | - | - | - | □ | - | - |
| Accessories | | | | | | | | |
| Firm Cable Installation | | | | | | | | |
| - Control Room | □ | □ | □ | □ | □ | □ | □ | - |
| - Treatment Room (RS232) | □ | □ | □ | □ | □ | □ | □ | - |
| Barometer/Thermometer | □ | □ | □ | □ | □ | □ | □ | □ |
| Computer Table/Trolley | □ | □ | □ | □ | □ | □ | □ | □ |

■ included □ optional - not available and/or not recommended

¹Part Nos for 230 V, Part Nos for 115 V: L981163, 64, 65 (BNT, TNC, M)



| PTW Detectors | Type No. | MP3 | MP3-M | MP3-T | MP3-XS SRS | MP3-XS IORT | MP3-P | MP2 | MP1 |
|--|--------------------|--|---|---|---|---|---|---|--------|
| | | ¹ L981160 BNT L981161 TNC L981162 M | L981166 BNT L981167 TNC L981168 M | L981199 BNT L981183 TNC L981182 M | L981172 BNT L981173 TNC L981174 M | L981383 BNT L981384 TNC L981385 M | L981185 BNT L981186 TNC L981233 M | L981175 BNT L981176 TNC L981177 M | T41019 |
| Semiflex Chamber 0.125 cm ³ | 31010 | ■ | ■ | - | ■ | ■ | - | ■ | □ |
| Semiflex Chamber 0.3 cm ³ | 31013 | □ | □ | - | □ | - | - | □ | □ |
| Farmer Chamber 0.6 cm ³ | 30013 | □ | □ | - | □ | □ | □ | □ | □ |
| Adv. Markus Chamber 0.02 cm ³ | 34045 | □ | □ | - | □ | □ | - | □ | □ |
| Markus Chamber 0.055 cm ³ | 23343 | □ | □ | - | □ | - | □ | □ | □ |
| Roos Chamber 0.35 cm ³ | 34001 | □ | □ | - | □ | ■ | □ | □ | □ |
| PinPoint 3D Chamber 0.016 cm ³ | 31016 | □ | □ | □ | □ | □ | - | □ | - |
| PinPoint Chamber 0.015 cm ³ | 31014 | □ | □ | - | □ | □ | - | □ | - |
| PinPoint Chamber 0.03 cm ³ | 31015 | □ | □ | ■ | □ | □ | - | □ | - |
| microLion Chamber 0.002 cm ³ | 31018 | □ | □ | - | □ | - | - | □ | - |
| Diamond Detector (1...6) mm ³ | 60003 ² | □ | □ | - | □ | - | - | □ | - |
| Dosimetry Diode P | 60016 | □ | □ | - | □ | □ | - | □ | - |
| Dosimetry Diode E | 60017 | □ | □ | - | □ | □ | - | □ | - |
| Dosimetry Diode SRS | 60018 ³ | □ | □ | - | □ | □ | - | □ | - |
| Bragg Peak Chamber 10.5 cm ³ | 34070 | - | - | - | - | - | □ | - | - |
| Thin-window Bragg Peak Chamber 10.5 cm ³ | 34080 | - | - | - | - | - | □ | - | - |
| Bragg Peak Chamber 2.5 cm ³ | 34073 | - | - | - | - | - | □ | - | - |
| SFD Mammo Chamber 6 cm ³ | 34069 | - | - | ■ | - | - | - | - | - |

■ included □ optional - not available and/or not recommended

¹Part Nos for 230 V, Part Nos for 115 V: L981163, 64, 65 (BNT, TNC, M); ²available for connecting system M only; ³small photon fields (10 cm x 10 cm) with max. 6 MV