



BEAMSCAN
Safety-Tag



◀ Send me back!

We will check for you whether anything has happened to your BEAMSCAN during delivery.

 PTW
THE
DOSIMETRY
COMPANY



1. Turn off the **BEAMSCAN Safety-Tag** by pressing and holding the center twice for three seconds each time. You will receive a confirmation through a green and red flashing of the LED.
 2. Remove the **BEAMSCAN Safety-Tag** from the BEAMSCAN by twisting the black screw on top.
 3. Contact us: shipping@ptwdosimetry.com
We will take care of the shipping of your **BEAMSCAN Safety-Tag**.
- Important: We need your collection address (as detailed as possible) and contact person as well as phone number and e-mail address.
4. Place the **BEAMSCAN Safety-Tag** into the enclosed package and attach the return label that you received from our shipping-team on top of the package.
 5. Send the package back to PTW Freiburg GmbH.